



**Room Setup and Standard Equipment:**

1. Standard Upper Endoscopy Setup
  - a. Endoscopy Tower (e.g., Olympus I40 Series or newer)
  - b. Gastroscope
  - c. Endoscope, 5 mm diameter or narrower and at least 105 cm in length (e.g., Olympus GIF-NI80 Gastroscope)
  - d. Room Suction
2. General Anesthesia Equipment
3. Laparoscopic Tower; incl. High Flow CO2 Insufflator

**General Supplies, Disposable:**

1. Standard drapes and patient accessories for upper G.I. endoscopy
2. ET Tube, IV Set for Anesthesia
3. Bite block, at least 20 mm internal diameter (e.g., Omni Bite Block: GI Supply GIS-11)
4. Astroglide and KY Lubricant
5. 20cc luer syringe, Emesis Basin, Gauze Sponges
6. Laparoscopic Insufflation Tubing, Standard

**Endoscopic Instruments, Reusable or Disposable:**

1. Straight scissors (e.g., Olympus FS-3L-1, Surgical Scissors 165cm x 2.8mm)
2. Biopsy forceps



USGI Medical  
1140 Calle Cordillera  
San Clemente, CA 92673

Customer Service  
1-866-788-USGI (8744)  
9:00 AM TO 5:00 PM Pacific Time